


## LACRIMAL DUCT RECONSTRUCTION INSTRUMENTS

Medtronic Solan offers a comprehensive selection of lightweight Intubation Instrument sets for reconstruction of the lacrimal outflow system, dacryocystorhinostomy, fistulization of the outflow system to the tear reservoir and other lacrimal system reconstructions.

### Fiberoptic Canalculus Intubation Set




85-90550 Two 15.5cm long, stainless steel and fiber probes, combined with provided penlight, provide direct visualization during intubation of the nasolacrimal system. 1 set per box.

## OPHTHO-BURR SYSTEM

Battery powered design aids in the removal of corneal foreign bodies and corneal rust rings. Uses 1.0mm, 2.5mm or 5.0mm burs.


- 85-50820 Complete system Includes one handpiece (85-50815), and one box of ten 1.0mm burs (85-50810).
- 85-50815 Handpiece only
- 85-50810 1.0mm bur, single use, sterile. Box of 10.
- 85-50850 2.5mm diamond bur, reusable.
- 85-50860 5.0mm diamond bur, reusable.

### Guibor Canalculus Intubation Set



85-90350 Two 15.5cm long, 0.051cm diameter, malleable, stainless steel probes, connected by 33cm of silicone tubing. Small handles on probes. 5 sets per box.

### C-Line Canalculus Intubation Set



85-9045 Two 17.7cm long, 0.064cm diameter, malleable, stainless steel probes connected by 29cm of silicone tubing. 5 sets per box.

### Ophtho-Burr System

